

Autograph AG-X Series

Shimadzu Precision Universal Tester



Autograph AG-X Series

Shimadzu Precision Universal Tester

Reliable, stress-free workflow

Reliability

Experience unsurpassed reliability and operability with the new Shimadzu Autograph AG-X series. This revolutionary new model delivers high-level controlled measurement performance utilizing a combination of patented technology from earlier models and world-class innovations. Development of this series focused on intuitive operation and convenient support functions, resulting in standard equipment that includes both a color TFT touch panel screen for PC-free operation and the very popular Smart Controller operation featured on many of our earlier testing machines. The TRAPEZIUM X software has also been newly developed, using cutting-edge Microsoft.Net technology to enhance user-friendliness. Shimadzu's new AG-X series advances testing in three areas performance, operability and support.

01 Superior Performance pg. 4

High reliability ensures complete data collection. Confidently perform comparisons with unknowns.

02 Unsurpassed Ease of Use pg. 6

Easy-to-use functions ensure smooth, trouble-free testing.

03 Quest for Convenience pg. 8

TRAPEZIUM X streamlines testing and eliminates confusion.

01 Superior Performance



High Control resolution results in testing reliability

A high-level control resolution of 0.0208 microns ensures consistent delivery of desired testing results.

Easy control of stress and strain

Auto tuning of control parameters is now possible in real time, based on measured test force and strain data. Comparisons can be safely made with unknown sample data, without the need for preliminary tests.

Accurate S-S Curves are achieved with highly precise load cells

Improve testing efficiency and ensure that virtually all of your testing can be performed without switching the load cell or jig, as a result of the wide, guaranteed load cell precision range of 1/1000 to 1/1.

Ultrahigh-speed sampling ensures no missed strength changes

Ultrahigh-speed 0.2 msec. (5 kHz) sampling ensures that sudden test force changes often seen at the start of testing can all be recorded. Easily change sampling condition settings during testing to investigate important regions in detail.

Convenient testing of actual objects

Up to 12 channels of data can be simultaneously read by a PC for immediate analysis. Test force readings along with data from multiple strain gauges may be collected during testing of actual objects.

Highly rigid frame guarantees Safety and reliability

The use of a crosshead guide provides high torsional rigidity to reduce movement and vibration of the machine, Shimadzu is accredited (JIS Q 17025, ISO/IEC 17025), based on JCSS as a calibration agency for uniaxial testing machines. JCSS cross-certifies with America's NVLAP and other certification standards throughout the world, providing world-wide traceability. (MRA certificate is offered by Shimadzu Shikenki Engineering Co., Ltd.)

· These products are CE compliant.

· AG-X units are manufactured by professionals at ISO9001 certified factories, ensuring years of worry-free operation.

02 Unsurpassed Ease of Use



6

Control at your fingertips - easy operation and data confirmation

With the Smart Controller, confirmation of test force and position real-time data is at your fingertips. Easily perform Start, Stop and other basic operations via this controller, and use the convenient jog wheel to adjust jig position in fine increments during bending and compression. You can even open or close the air chucks during tensile tests and operate the automatic extensometer.



An optional LCD touch panel means you can guickly select testing methods without having to connect a PC. Easily view graphs of data directly on the LCD screen.

Store test methods in USB memory

After storing testing methods on a USB memory device, simply insert the device into the testing machine to perform testing without a PC. Measurement data can also be automatically saved to USB memory. After testing, bring your USB memory device back to your office PC to analyze data and create reports. (Requires LCD touch panel and TRAPEZIUM X software.)

Self-diagnostics help cover all bases

Self-check function (12 items, including motor pulse, sensor amplifier, and board power supply) confirms that the instrument is in perfect working order. If desired, notification of pre-set maintenance periods is also possible. (Some check items require a special jig.)

Simple load cell installation (option)

Use this load cell quick attach/release unit with the table-top type 10 kN capacity AG-X unit, which usually requires frequent load cell changes.

AG-X units with a capacity of 20 kN or more can attach the optional small-capacity load cell attachment plate to the bottom of the crosshead, eliminating the need to detach the original load cell.



Load cell attach/release (10 kN type)



Floor model (with small-capacity load cell attached)

Safety equipment

· Safety cover

This cover is designed to control scattering of the test specimen during testing and the interlock improves operation safety.

Safety functions

If force changes exceed a certain level during specimen setting or return, the testing machine is stopped by the safety function.

· Dual emergency-stop switches

As a safety measure, emergency-stop switches are provided on both sides.



Safety cover (option) attached





03 Quest for Convenience





03 Quest for Convenience

Receive data quickly

1. Speed, dimension, and report information can be entered quickly and directly from the main window using the [Quick Panel].

2. Advanced navigation system with learning functions

• AG-X is equipped with a Navigation Bar that shows only the functions required for a selected situation. This allows you to efficiently perform continuous testing using simple, straightforward procedures and by pressing large, easy-to-read buttons.

AG-X is also equipped with a "Learning function" that records user actions for each situation and adds frequently-used functions to the Navigation Bar. This means that the more you use the machine, the better the "fit" is to your unique operation style, effectively speeding up your workflow.

- 3. Functions include re-test, file synthesis, as well as specimen insertion, addition and order changes in any position.
- Re-test: A portion of a batch test can be retested, and the prior test results replaced.
- Extra lot tests: batches (lots) can be added, increasing the total number of tests.
- A variety of setting changes are possible before and after testing. Specimens can be inserted in any position or added to only a specific batch, and the specimen order can be changed after testing is completed.







Generate detailed reports

Richly expressive report creation includes free positioning of report elements and a wealth of web-compatible output functions.



- **Reports can be output in PDF, Microsoft Word, Excel and HTML formats.** Output reports created with Report Designer in a wide variety of useful formats. (Charts and tables with ruled lines cannot be output in Word and HTML.) After export, use your everyday software to customize the report.
- WebPlus function (option)
 Installing the WebPlus option on your server PC allows reanalysis and printing via Internet Explorer, even on a PC not equipped with TRAPEZIUM X.



Choose from four software components to fit your specific application

TRAPEZIUM X includes four software components - Single, Cycle, Control and Texture. This allows you to purchase only the components that meet your specific testing needs. When multiple software components are purchased, easily switch between modes at a single touch, without starting up separate software.

fest Mode	_
Single	~
Single	2
Cycle	N
Control	5
Texture	w

Single software

Performs general single-direction testing. Examples include tensile, compression, bending and peeling tests.



Control software

Create any testing machine operation pattern. Perform foam rubber compression and holding cycle tests.



Cycle software

Similar to endurance testing, this software is used for testing where force is repeatedly applied and then released.



Texture software

Measures the features (texture) of foods and pharmaceuticals. Produce special data processing results, including mastication, jelly strength and adhesion.



Accessories

[Experience the range of possibilities available with this full-featured system]

Accessories Lineup

Tensile tests

Combine grips and extensometers with the testing machine.

Grips

Used to grip the sample, a wide variety is available to accommodate different specimen types and test force amounts.

Non-shift wedge type grips <MWG> Applications: Plastics, Metals, Wood

Grip		Standard grip face			Upper grip capacity		
capacity	Grip face	Clearance (mm)	Grip width (mm)	Grip length (mm)		Part No.	
250kN		0 to 8.5	50	75	33	343-07979-12	
100kN	File teeth for	0 to 7	40	55	10	346-52791-12	
50kN		0 to 7	40	55	9.5	346-52791-11	
20kN	flat specimens	0 to 7	25	55	3.6	346-52653-12	
5kN		0 to 7	25	55	3.6	346-52653-11	

Screw type flat grips <SCG> Applications: Rubber, Plastics, Textiles, Cloth, Paper

Grip		Standard grip face			Upper grip capacity	
capacity	Grip face	Clearance (mm)	Grip width (mm)	Grip length (mm)	(kg)	Part No.
5kN	File teeth	0 to 16	60	50	2	345-52326-04
1kN	File leelii	0 to 15	50	30	0.7	346-52327-04
50N	Flat	0 to 14	35	25	0.3	346-52328-04

Pneumatic flat grips <PFG> Applications: Rubber, Plastics, Textiles, Cloth, Paper

Grip	External din	nensions (mm)	Crin width (mm)	Clearance (mm)	Upper grip capacity	
capacity	w	L (upper/lower)	anp watn (mm)	Clearance (mm)	(kg)	Kit No.*1, *2
10kN	154	268.5 / 278.5	60	0 to 10	<u> </u>	346-53916-XX
5kN	154	224 / 235	60	0 to 6	5.7	346-53849-XX
1kN	102	163 / 174	50	0 to 6	1.7	364-53848-XX
50N	64	118 / 135	35	0 to 6	0.4	346-53847-XX

*1 Grips with foot-valve units and crosshead-linked control functions are also available.

*2 Grips can be opened and closed via the Smart Controller when using the crosshead-linked control kit.

Extensometers

Extensometers improve elongation measurement accuracy.

Strain gauge type one-touch extensometer <SSG-H Series>

SSG-H series extensometers conform to JIS B7741 Class 0.5 and JIS K7161 (SSG 50-10SH only). They can be attached using just one touch.

Model	Gauge length (mm)	Measuring range (mm)	Kit No.
SSG25-50H	25	12.5 5.25 2.5 1.24	346-53875-23
SSG25-100H	25	25 12.5 5 2.5	346-53875-24
SSG50-10H	50	5 2.5 1 0.5	346-53875-51
SSG50-10SH	50	5 2.5 1 0.5	346-53875-56

* Calibration cables (for SGI) are included with each kit.

* Precision is JIS B7741 Class 0.5 or Class 1, depending on the conditions.

DVE series non-contact video extensometer

DVE series extensometers use two cameras to provide a wide measurement range and high measurement precision.

- GL50 mm Maximum measurement range 30 mm Absolute precision of indicated value $\pm 3 \ \mu$ m (JIS B7741 Class 1) Relative precision $\pm 1\%$
- GL20 mm Maximum measurement range 140 mm Absolute precision of indicated value $\pm 6 \ \mu$ m (JIS B7741 Class 2) Relative precision $\pm 1\%$

(Elongation is measured via the PC monitor screen. Elongation measurement in an environmentally-sealed chamber is also possible.)





Screw type flat grips







Compression tests

Simply attach the compression plate kit to the main unit to perform compression testing.

· Compression plate kit Applications: Plastics, Metals, Rubber, Wood, Cement

Used to compress the specimen, several types are available to accommodate different specimen types and test force amounts.

Fixed type

Kit No.	Maximum capacity	Upper plate dimensions (mm) diameter by thickness	Upper plate mass (kg)	Operational temperature (°C)
346-53882-XX		ø100 x 25	1.6	K
346-53884-XX	250kN	ø50 x 25	0.5	0 to 40
346-53885-XX		ø200 x 40	6.3	

Spherical seat type

Kit No.	Maximum capacity	Upper plate dimensions (mm)	Upper plate mass (kg)	Operational temperature (°C)
346-53883-XX	250kN	ø100	3.8	0 to 40

* With spherical compression plates, only the upper plate is spherical.

Spherical seat-type compression plates provide contact flexibility for uniform load application.

* Select the kit number that corresponds to the load cell used.



Fixed type compression plates



Spherical seat type compression plates

Bending tests Simply attach the bending test jig kit to the main unit to perform bending testing.

Bending test kit

Select the kit number appropriate for the load cell used.

Kit No.	Max. test force	Punch tip radius x width (mm)	Support tip radius x width (mm)	Support spacing (mm)	Operational temperature (°C)	Applicable test standards
346-53887-XX	10kN	R5 x 34	R2 x 34	20 to 200	lan.	JIS K6911, JIS K6902*1, JIS C6481*2, JIS K7171, ISO 178, Specimens with thickness of 3 mm or less
340-33007-77	TOKIN	****	R5 x 34	2760	Xa.	JIS K7171, ISO 178, Specimens with thickness above 3 mm
F CD A C		R1/8" x 72	R1/8" x 110	0.8 to 8"	0 to 40	ASTM D790 (Test method 1)
		R5 x 72	R2 x 110	50 to 500	EXE	JIS K6911, JIS K6902*1, JIS C6481*2, JIS K7171, ISO 178, Specimens with thickness of 3 mm or less
346-53888-XX	100kN		R5 x 110		THE PAR	JIS K7171, ISO 178, Specimens with thickness above 3 mn
	100 V	R1/8" x 72	R1/8" x 110	2 to 20"	VALUE V	ASTM D790 (Test method 1*3)

*1 Corresponds to bending strength. Compatible with support spacing from 20 mm to 200 mm.

*2 Corresponds to bending strength. *3 Compatible with support spacing from 2 inches to 20 inches.

When the SIE or SES extensometer is used, the following adaptor is required. 346-55658-XX

Adhesion test

• Adhesive tape peeling test device Applications: Plastics, Rubber

Specimen table slides in accordance with upper grip movement to maintain a 90 degree peeling angle.

Upper grip: 1 kN Flat screw type, 1 unit

Kit No.	Capacity	Applicable specimen (width x thickness mm)	Operational temperature (°C)	Applicable test standards
346-53865-XX	1kN	50 x 5 to 2	-10 to +60	JIS Z0237 Adhesive tape Adhesive tape test method (90° peeling test) JIS Z1528 Double-sided adhesive tape adhesion

Mechanism differs from JIS example.

Select the kit number appropriate for the load cell used.

For details on test jigs not listed in this catalog, please refer to the separate Accessories catalog.



3-point bending test of plastic specimen

Accessories

Accessories Lineup

- Various grips
- Adhesion test devices
- Shearing test devices
- Tear test devices
- Devices for needle insertion resistance
- measurement
- Flow test devices

- Friction coefficient measuring devices
- Plastic bearing strength test devices
- Deep-drawing test devices
- · Nail withdrawal resistance test devices (from wood)
- Lumber hardness test devices
- · Lumber cleavage test devices

- · Powder molding properties test devices
- · Controlled atmosphere test devices
- Displacement measuring devices Test force measuring devices Others

Grips and devices for testing actual objects

 Pneumatic automatic grips **PWG** series

Pneumatic capstan type grips

force can be maintained.

Specimens such as threads and cords

are held by the capstan. Initial tensile

Air motor opens and closes grips to shorten time required for testing.



- Spring tensile test jig
- Efficiently evaluates mechanical characteristics of tensile coil springs.



• Printed Circuit Board 45-degree peeling test jig Used for peeling tests of electronic



parts on a printed circuit board.

Auto Extensometer

SIE-560S

This extensometer uses a highprecision strain-gauge sensor and magnetic induction sensor to automatically set the gauge marker positions. The extensometer can be automatically attached/removed.

When used with bending test jigs, the following adaptor is required. 346-55658-XX

Soft material extensometer **SES-1000**

Easily and accurately measures large elongation amounts. (PAT No. 3724136)

When used with bending test jigs, the following adaptor is required. 346-55658-XX



 Compression plate displacement measurement device Measures displacement of compression plates during compression tests.



- Strain gauge type width sensor Measures changes in specimen width.





 Jigs for cyclic bending tests of printed circuit boards

Allows cyclic bending testing of printed circuit boards containing parts. Combining this jig with the optional resistance meter increases test efficiency by quickly detecting internal damage and stopping the test.

- Silicon chip (die) 3-point bending Max. capacity: 500 N

Punch dimensions: tip R0.3 x 20 mm Support dimensions: tip R0.3 x 20 mm Span: 1 mm to 20 mm

Conforms to SEMI G86-0303 * SEMI: Semiconductor Equipment and Materials International

Controlled atmosphere test device

 Thermostatic chamber **TCE** series

test jig

This compact chamber enables testing across a wide temperature range of -70 °C to +280 °C.

 Bellows-type long stroke thermostatic chamber Bellows design is well-suited to testing of highly elastic materials.





- Long-span 4-point bending test device for wooden specimens Capacity: 100 kN Punch: width 1320 mm Punch distance: 100 mm to 1350 mm Support: width 1320 mm Support distance: 100 mm to 4000 mm
- PC card insertion and removal test jig Allows repeated insertion and removal of PC cards (memory cards) and PC card adaptors.





Device for torsion test in thermostatic chamber Temperature range: -60 °C to +250 °C Torsion capacity: 300 N-m Torsion speed: 1 rpm to 0.01 rev/min



 In-chamber tensile test device This specially-designed device passes the crosshead through a thermostatic chamber to assure a long effective stroke.



Specifications

[Table-Top AG-X]

1. Model Name		AG-10N / 20NX AG-50N / 100NX	AG-500N / 1kNX	AG-5kNX	AG-10kNX	AG-20kN / 50kNXD		
2. Capacity		10N / 20N 50N / 100N	500N / 1kN	5kN	10kN	20kN / 50kN		
. Loading Metho	od VA	Direct, hig	h-precision, constant-rate	strain control using non-ba	acklash precision ball-screv	w drive		
Pre	High-precision unit 1/1000	Not available Within ± 0.3% Conforms to	of displayed test force (for	or 1/100 to 1/1000 of load o or 1/1 to 1/100 of load cell 10002-2 Grade 0.5, ISO 75	cell rated capacity) rated capacity) 500-1 Class 0.5, BS1610 Cl	ass 0.5,		
4. Test Force Measurement Standard-precision unit		Within ± 1% of displayed test force (for 1/1 to 1/1000 of the load cell rated capacity) Conforms to JIS B7721 Class 1, EN 10002-2 Grade 1, ISO 7500-1 Class 1, BS1610 Class 1, DIN51221 Class 1, and ASTM E4*3						
	1/500	Within ± 1% of displayed test force (for 1/1 to 1/500 of load cell rated capacity) Conforms to JIS B7721 Class 1, EN 10002-2 Grade 1, ISO 7500-1 Class 1, BS1610 Class 1, DIN51221 Class 1, and ASTM E4*3						
	Test force calibration			nd compression forces calibra calibration of tensile force, cor	tion mpression force, or both tensil	e and compression for		
. Crosshead Sp	eed Bange			Free step-less setting				
mm/min)				0.0005 ~ 1000				
	Maximum Return Speed		1:	500		1200		
. Crosshead Sp	eed Precision*1			±0.1%				
. Position Contr	rol Resolution		0.0	25µm		0.0208µm		
3. Crosshead Spee	ed and Allowed Test Force		Maxim	um load capacity for all spe	eeds			
. Crosshead-Ta Tensile stroke) '	ble Clearance (mm) *2	Max. 1150 mm (850 mm)	Max. 1150 mm (850 mm)	Max. 1150 mm (780 mm)	Max. 1150 mm (600 mm)	Max. 1060 mm (655 mm): 20 kl (605 mm): 50 kl		
0. Effective Tes	t Width (mm)		4	20		500		
		1.4.1.1.4.				-		
1. Crosshead Population				coder measurement, digita				
	display methods	Within :			I display indicated value is below 10	mm		
Detection	osition display methods Precision	Within :				mm		
Detection 2. Data Capture	osition display methods Precision Rate	Within :		nowever, ±0.01 mm when i		mm		
Detection 12. Data Capture 13. Data Samplin	osition display methods Precision Rate Ig Rate	Within :	£0.1% of indicated value, I	nowever, ±0.01 mm when i		mm 120		
11. Crosshead Pr Detection 12. Data Capture 13. Data Samplin 14. Frame Rigidit	osition display methods Precision Ing Rate ty (kN/mm)	Within : • Automatic reading of load • Fine adjustment of crosshu • Test force and stroke disp • External analog output (2 ch • External digital input (2 ch • External digital input (2 ch • Internal amps - 4 ports (one is used for test force and • USB interface (for PC) / H • Recorder output (optional) • Dataletty output (optional) • Pneumatic grip interlock o	e0.1% of indicated value, I cell properties ead position ay channels) nannels) annels) annels) another for analog input) post interface (for USB memo	nowever, ±0.01 mm when i 5000Hz 300kHz 42 • Automatic test for • Test force auto a • Break detection • Crosshead spece • Stress value dis • Soft limit detection If only optional L • Single testing co Control of testing • PEAK and BRE/ • Method internal	indicated value is below 10 prce and strain control (with a zero / auto calibration	120 auto tuning) iisplay splay ad speed pre-setting		
Detection 2. Data Capture 3. Data Samplin 4. Frame Rigidit 5. Standard Fur	osition display methods Precision Precision Precision pg Rate ty (kN/mm)	Automatic reading of load Fine adjustment of crosshu Test force and stroke disp External analog output (2 cf External digital input (2 ch Internal amps - 4 ports (one is used for test force and USB interface (for PC) / H Recorder output (optional) Dataletty output (optional)	cell properties pad position ay phannels) annels) another for analog input) post interface (for USB memory post interface (for USB memory peration (optional)	nowever, ±0.01 mm when i 5000Hz 300kHz 42 • Automatic test for • Test force auto a • Break detection • Crosshead spece • Stress value dis • Soft limit detection If only optional L • Single testing co Control of testing • PEAK and BRE/ • Method internal	indicated value is below 10 proce and strain control (with a zero / auto calibration / auto return ad free setting / cycle count d play / extensometer value dis on / self diagnostics .CD touch panel is used: .CD touch pan	120 auto tuning) iisplay splay ad speed pre-setting		
Detection 2. Data Capture 3. Data Samplin 4. Frame Rigidit 5. Standard Fur	display methods osition Precision Image: Rate Ing Rate ty (kN/mm)	Automatic reading of load Fine adjustment of crosshe Test force and stroke disp External analog output (2 cf External digital input (2 ch Internal amps - 4 ports (one is used for test force and USB interface (for PC) / H Recorder output (optional) Dataletty output (optional) Pneumatic grip interlock o	cell properties pad position ay channels) annels) another for analog input) post interface (for USB memory peration (optional) For 500N/1 kN	however, ±0.01 mm when i 5000Hz 300kHz 42 • Automatic test fo • Test force auto z • Break detection • Crosshead spee • Stress value dis • Soft limit detecti If only optional L • Single testing cc cory) Control of testing • PEAK and BRE/ • Method internal • Japanese/Englis	indicated value is below 10 proce and strain control (with a zero / auto calibration / auto return ad free setting / cycle count d play / extensometer value dis on / self diagnostics .CD touch panel is used: .CD touch panel is used: .Dtrol / Cycle testing control / g conforming to standards AK values display / Crosshea memory file (20 files) sh switchover / S-S curve dis	120 auto tuning) isplay splay ad speed pre-setting		
Detection 2. Data Capture 3. Data Samplin 4. Frame Rigidit 5. Standard Fur 6. Accessories	display methods osition Precision ing Rate ing Rate ty (kN/mm)	Automatic reading of load Fine adjustment of crosshe Test force and stroke disp External analog output (2 cf External digital input (2 ch Internal amps - 4 ports (one is used for test force and USB interface (for PC) / H Recorder output (optional) Dataletty output (optional) Pneumatic grip interlock o For 10 N/20 N For 50 N/100 N	E0.1% of indicated value, I cell properties ead position ay channels) annels) another for analog input) ost interface (for USB memory peration (optional) For 500N/ 1 kN For 500N/ 1 kN	nowever, ±0.01 mm when i 5000Hz 300kHz 42 • Automatic test for • Test force auto : • Break detection • Crosshead spee • Stress value dis • Soft limit detecti If only optional L • Single testing cc ory) • Control of testing • PEAK and BRE/ • Method internal • Japanese/Englis For 5 kN For 5 kN	indicated value is below 10 proce and strain control (with a zero / auto calibration / auto return ed free setting / cycle count d play / extensometer value dis on / self diagnostics .CD touch panel is used: .CD touch	120 auto tuning) isplay splay ad speed pre-setting play For 20 kN/50 kl		
Detection 2. Data Capture 3. Data Samplin 4. Frame Rigidit 5. Standard Fur 6. Accessories	display methods osition Precision Precision Ing Rate ty (kN/mm) nctions Load cell CAL. cable Others	Automatic reading of load Fine adjustment of crosshe Test force and stroke disp External analog output (2 cf External digital input (2 ch Internal amps - 4 ports (one is used for test force and USB interface (for PC) / H Recorder output (optional) Dataletty output (optional) Pneumatic grip interlock o For 10 N/20 N For 50 N/100 N	E0.1% of indicated value, I cell properties ead position ay channels) annels) another for analog input) ost interface (for USB memory peration (optional) For 500N/ 1 kN For 500N/ 1 kN Standard tool s	nowever, ±0.01 mm when i 5000Hz 300kHz 42 • Automatic test for • Test force auto 2 • Break detection • Crosshead spee • Stress value dis • Soft limit detecti If only optional L • Single testing co cory) • PEAK and BRE/ • Method internal • Japanese/Englis For 5 kN Et, instruction manuals, lim	indicated value is below 10 proce and strain control (with a zero / auto calibration / auto return ed free setting / cycle count d play / extensometer value dis on / self diagnostics .CD touch panel is used: .CD touch	120 auto tuning) isplay splay ad speed pre-setting play For 20 kN/50 kI For 20 kN/50 kI		
Detection 2. Data Capture 3. Data Samplin 4. Frame Rigidit 5. Standard Fur 6. Accessories 7. Dimensions	display methods osition Precision Precision Ing Rate ty (kN/mm) nctions Load cell CAL. cable Others	Automatic reading of load Fine adjustment of crosshe Test force and stroke disp External analog output (2 cf External digital input (2 ch Internal amps - 4 ports (one is used for test force and USB interface (for PC) / H Recorder output (optional) Dataletty output (optional) Pneumatic grip interlock o For 10 N/20 N For 50 N/100 N	E0.1% of indicated value, I cell properties ead position ay channels) annels) another for analog input) ost interface (for USB memory peration (optional) For 500N/ 1 kN For 500N/ 1 kN Standard tool s	nowever, ±0.01 mm when i 5000Hz 300kHz 42 • Automatic test for • Test force auto : • Break detection • Crosshead spee • Stress value dis • Soft limit detecti If only optional L • Single testing cc ory) • Control of testing • PEAK and BRE/ • Method internal • Japanese/Englis For 5 kN For 5 kN	indicated value is below 10 proce and strain control (with a zero / auto calibration / auto return ed free setting / cycle count d play / extensometer value dis on / self diagnostics .CD touch panel is used: .CD touch	120 auto tuning) isplay splay ad speed pre-setting play For 20 kN/50 kl		

*1: Crosshead speed precision is calculated using crosshead transfer amount within a specified period of time for the crosshead speed of 0.5 mm/min to 500 mm/min under normal conditions.

*2: Tensile stroke is the value used when attaching the MWG (non-shift wedge type) grips. Stroke can be extended.

Values under 5 kN are with SCG (screw type flat) grips attached. *3:JIS B7721, EN 10002-2, ISO 7500-1, and ASTM E4 standards recommend re-verification after installation of testing machine.

*4:In conformity with CE Mark regulation

*5:Displayed units are by default expressed in SI system. Other unit system is selectable as Metric or English imperial by a keystroke.

*6:Software is available in several languages (English,Spanish,Chinese,Japanese etc.)

* Values in this catalog have been measured based on separately-approved test standards.



Up to 10 kN table-top model

· Laptop PC and table are optional.

• Table used in image (for up to 10 kN table-top model) is not a Shimadzu product.

50 kN table-top model

- · LCD touch panel, laptop PC, and table are optional.
- Table used in image is not a Shimadzu product.

Installation Space

(Dimensions given for left, right and back of main unit are space required for maintenance.)



Model	Mass (approx. kg)	Power Requirement - consumed power is in ()	Installation Environment
AG-10NX to 10kNX	135	Single phase 100 to 100/115 to 130/ 220 to 230/240 V (switching system) 50 to 60 Hz 1.5 kVA (450 W)	Temp.: 5 °C to 40 °C Humidity: 20% to 80% (no condensation) Voltage fluctuation: ±10% max. - Vibration: Frequency 10 Hz max. Amplitude 5 µm max.
AG-20 / 50kNX	162	Single phase 200 to 230 V 50 to 60 Hz 5 kVA (1.2 kW)	

(NB) Grounding of 100W or more is required.

AG-X Capacity/Model Types and Kit Numbers

Kit Number (AG-X unit + load cell set + upper/lower joints)

346-567XX- X1 Capacity Testing force measurement rating and guaranteed range 05 : 50kN table-top model 06 : 20kN table-top model 5 : Class 1 1/500 6 : Class 1 1/1000 07:10kN 08:5kN 7 : Class 0.5 1/1000 09 1kN 10 : 500N 11 : 100N 12: 50N 13: 20N (Ex.) Capacity 5 kN, test force measurement rating Class 1 14: 10N

LCD touch panel unit (for table-top models) P/N: 346-55227-51

Kit Number (AG-X unit + load cell set + upper/lower joints)

Load cell one-touch attachment unit (for 10 kN or less table-top models)

P/N: Load cell one-touch attachment unit 346-55042 346-55042-01 Load cell attachment

The Load cell one-touch attachment unit includes a load cell attachment. One load cell attachment is required for each load cell attached.

Specifications

[Floor Type AG-X]

1. Model Name			Floor Type AG-X					
n model Name		AG-20kN / 50kNXD	AG-100kNX	AG-250kN / 300kNX				
2. Capacity		20kN / 50kN	100kN	250kN / 300kN				
3. Loading Meth	od YA	Direct, high-precision, constan	Direct, high-precision, constant-rate strain control using non-backlash precision ball-screw drive					
Pre	High-precision unit 1/1000 (1/250 for 250 kN and 300 kN models)	Within ± 0.5% of displayed test force (for 1/100 to 1/1000 of load cell rated capacity) Within ± 0.3% of displayed test force (for 1/1 to 1/100 of load cell rated capacity) Conforms to JIS B7721 Class 0.5, EN 10002-2 Grade 0.5, ISO 7500-1 Class 0.5, BS1610 Class 0.5, DIN51221 Class 1, and ASTM E4*3Within ± 0.5% of displayed test for 						
Precision 4. Test Force Measurement Standard-precision unit		Within ± 1% of displayed test force (for 1/1 to 1/1000 of load cell rated capacity) Conforms to JIS B7721 Class 1, JIS B7733 Class 1, EN 10002-2 Grade 1, ISO 7500-1 Class 1, BS1610 Class 1, DIN51221 Class 1, and ASTM E4*4						
Measurement	1/500	Within ± 1% of displayed test force (for 1/1 to 1/500 of the load cell rated capacity) Conforms to JIS B7721 Class 1, EN 10002-2 Grade 1, ISO 7500-1 Class 1, BS1610 Class 1, DIN51221 Class 1, and ASTM E4*3						
	Test force calibration	Automatic calibration Standard-precision type: Tens High-precision type: Choose		on force, or both tensile and compression forces				
5. Crosshead Sp	and Pango		Free step-less setting					
(mm/min)		0.0005 to 10	000	0.0005 to 500				
	Maximum Return Speed	1200						
6. Crosshead Sp	eed Precision*1		±0.1%					
7. Position Cont	rol Resolution	0.0208µn	1	0.0104µm				
8. Crosshead Spe	ed and Allowed Test Force	Maximum load capacity	r for all speeds	0.0005 ~ 250mm/min : 300kN 0.0005 ~ 500mm/min : 250kN				
9. Crosshead-Ta (Tensile stroke)	ble Clearance (mm) *2	Max. 1265 mm (850 mm): 20 kN (800 mm): 50 kN	Max. 1250 mm (760 mm)	Max. 1440 mm (600 mm)				
10. Effective Tes	st Width (mm)		595					
11. Crosshead P Detection	Measurement and display methods osition Precision							
	riccision	Within ±0.1% of indicated value, but ±0.01 mm when the indicated value is below 10 mm						
12. Data Capture	Rate	5000Hz						
13. Data Samplin	ng Rate	1124112	300kHz					
14. Frame Rigidi								
	ty (kN/mm)	175	300	400				
	ty (kN/mm)	 Automatic reading of load cell properties Fine adjustment of crosshead position Test force and stroke display External analog output (2 channels) 	 Automatic test force Test force auto zero Break detection / au Crosshead speed fr 	and strain control (with auto tuning) / auto calibration to return ee setting / cycle count display				
15. Standard Fur		 Automatic reading of load cell properties Fine adjustment of crosshead position Test force and stroke display 	Automatic test force Test force auto zerc Break detection / au Crosshead speed fr Stress value display Soft limit detection / If only optional LCD Single testing control PEAK and BREAK v Method internal mer	and strain control (with auto tuning) / auto calibration to return ee setting / cycle count display / extensometer value display self diagnostics touch panel is used: ol / Cycle testing control / informing to standards /alues display / Crosshead speed pre-setting				
15. Standard Fu		 Automatic reading of load cell properties Fine adjustment of crosshead position Test force and stroke display External analog output (2 channels) External analog input (2 channels) External digital input (2 channels) Internal amps - 4 ports (one is used for test force and another for analog input USB interface (for PC) / Host interface (for USB Recorder output (optional) Dataletty output (optional) 	Automatic test force Test force auto zerc Break detection / au Crosshead speed fr Stress value display Soft limit detection / If only optional LCD Single testing control PEAK and BREAK v Method internal mer	and strain control (with auto tuning) / auto calibration to return ee setting / cycle count display / extensometer value display self diagnostics touch panel is used: o/ / Cycle testing control / informing to standards /alues display / Crosshead speed pre-setting mory file (20 files)				
	nctions	 Automatic reading of load cell properties Fine adjustment of crosshead position Test force and stroke display External analog output (2 channels) External analog input (2 channels) External digital input (2 channels) Internal amps - 4 ports (one is used for test force and another for analog input USB interface (for PC) / Host interface (for USB Recorder output (optional) Dataletty output (optional) Pneumatic grip interlock operation (optional) 	 Automatic test force Test force auto zerc Break detection / au Crosshead speed fr Stress value display Soft limit detection / If only optional LCD Single testing control Control of testing control of testing	and strain control (with auto tuning) / auto calibration to return ee setting / cycle count display / extensometer value display self diagnostics touch panel is used: ol / Cycle testing control / onforming to standards /alues display / Crosshead speed pre-setting mory file (20 files) witchover / S-S curve display				
	nctions	Automatic reading of load cell properties Fine adjustment of crosshead position Test force and stroke display External analog output (2 channels) External analog input (2 channels) External digital input (2 channels) Internal amps - 4 ports (one is used for test force and another for analog input USB interface (for PC) / Host interface (for USB Recorder output (optional) Dataletty output (optional) For 20 kN/50 kNX For 20 kN/50 kNX	Automatic test force Test force auto zerc Break detection / au Crosshead speed fr Stress value display Soft limit detection / If only optional LCD Single testing control Control of testing co PEAK and BREAK Method internal mer Japanese/English s For 100 kNX	and strain control (with auto tuning) and strain control / display and strain control / display / Crosshead speed pre-setting mory file (20 files) witchover / S-S curve display For 250 kN/300 kN For 250 kN/300 kN				
16. Accessories	Load cell CAL. cable Others	Automatic reading of load cell properties Fine adjustment of crosshead position Test force and stroke display External analog output (2 channels) External analog input (2 channels) External digital input (2 channels) Internal amps - 4 ports (one is used for test force and another for analog input USB interface (for PC) / Host interface (for USB Recorder output (optional) Dataletty output (optional) For 20 kN/50 kNX For 20 kN/50 kNX	Automatic test force Test force auto zerc Break detection / au Crosshead speed fr Stress value display Soft limit detection / If only optional LCD Single testing control Control of testing cc PEAK and BREAK Method internal mer Japanese/English s For 100 kNX For 100 kNX	and strain control (with auto tuning) and strain control display are setting / cycle count display self diagnostics touch panel is used: bl / Cycle testing control / unforming to standards values display / Crosshead speed pre-setting mory file (20 files) witchover / S-S curve display For 250 kN/300 kN For 250 kN/300 kN ning label				
	Load cell CAL. cable Others	Automatic reading of load cell properties Fine adjustment of crosshead position Test force and stroke display External analog output (2 channels) External analog input (2 channels) External digital input (2 channels) Internal amps - 4 ports (one is used for test force and another for analog input USB interface (for PC) / Host interface (for USB Recorder output (optional) Dataletty output (optional) For 20 kN/50 kNX For 20 kN/50 kNX	Automatic test force Test force auto zerc Break detection / au Crosshead speed fr Stress value display Soft limit detection / If only optional LCD Single testing control Control of testing cc PEAK and BREAK Method internal mer Japanese/English s For 100 kNX For 100 kNX tool set, instruction manuals, limit warr	and strain control (with auto tuning) and strain control display are setting / cycle count display self diagnostics touch panel is used: bl / Cycle testing control / anforming to standards values display / Crosshead speed pre-setting mory file (20 files) witchover / S-S curve display For 250 kN/300 kN For 250 kN/300 kN				

*1: Crosshead speed precision is calculated using crosshead transfer amount within a specified period of time for the crosshead speed of 0.5 mm/min to 500 mm/min under normal conditions.

*2: Tensile stroke is the value used when attaching the MWG (non-shift wedge type) grips. Stroke can be extended.

Values under 5 kN are with SCG (screw type flat) grips attached.

*3:JIS B7721, EN 10002-2, ISO 7500-1, and ASTM E4 standards recommend re-verification after installation of testing machine.

*4: In conformity with CE Mark regulation

*5:Displayed units are by default expressed in SI system.

Other unit system is selectable as Metric or English imperial by a keystroke. *6:Software is available in several languages (English,Spanish,Chinese,Japanese etc.)

* Values in this catalog have been measured based on separately-approved test standards.



20 kN, 50 kN floor model

· Laptop PC and table are optional.

Installation Space

(Dimensions given for left, right and back of main unit are space required for maintenance.)



100 kN floor model

• LCD touch panel, laptop PC, and table are optional.

AG-X Capacity/Model Types and Kit Numbers

Kit Number (AG-X unit + load cell set + upper/lower joints) 346-567XX-X1

Capacity	Testing force measurement rating and guaranteed range
00 : 300kN	5 : Class 1 1/500
01 : 250kN	6 : Class 1 1/1000
02 : 100kN	7 : Class 0.5 1/1000
03 : 50kN 04 : 20kN	(Ex.) Capacity 100 kN, test force measurement rating Class 1
	Guaranteed range of 1/500: 346-56702-51

LCD touch panel unit (for floor models) P/N: 346-55227-52

Model	Mass (approx. kg)	Power requirement - consumed power is in ()	Installation Environment
AG-20kN / 50kNX	620	Three phase 200 to 230 V 50 to 60 Hz 5 kVA (1.2 kW)	
AG-100kNX	800	Three phase 200 to 230 V 50 to 60 Hz 7 kVA (2.0 kW)	Humidity: 20% to 80% (no condensation)
AG-250kN / 300kNX	920	Three phase 200 to 230 V 50 to 60 Hz 7.5 kVA (2.5 kW)	Vibration: Frequency 10 Hz max. Amplitude 5 μ m max.

(NB) Grounding of 100W or more is required.

Extensions to the main unit [Table-top models and floor models]

Unit: mm

Wide frame series

Models with wider effective test widths (975 mm, 1100 mm and 1375 mm) than the standard type (595 mm) are also available for testing large-size actual object specimens. (Floor models only.)

Reinforced yoke series

Use this series when conducting tests between the crosshead and yoke.

Ultralow-speed crosshead model

The crosshead speed range can be widened to include extremely low speeds.

	Ultralow-speed models	
Speed Range (mm/min)	0.00005 mm/min to 1000 mm/min (250 kN and 300 kN models are limited to 500 mm/min.)	
Part Number	345-50522	

Large capacity series

With maximum capacities of 500 kN, 600 kN and 1000 kN, these models are used for large-capacity testing in heavy-industry fields such as steel, construction, and shipbuilding. They can be customized upon request.

Extended column models

Models with extended columns are useful for testing materials requiring long tensile strokes. (Clearance between the yoke and table is extended 250 mm, 500 mm, or 750 mm.)

High-speed return models

Table-top models less than 2 kN can be equipped with drives that have a return speed of 3000 mm/min, and a crosshead speed range of 0.001 mm/min to 2500 mm/min.

Testing and Evaluation Machines



(Shanghai) Co., Ltd.

Pte Ltd.

Instruments (Oceania) Pty Ltd.

Comercio Ltda.

*Windows® is a registered trademark of the Microsoft Corporation. *TRAPEZIUMX are designed for Windows® Me/2000/XP. They cannot be used with Windows® 3.1 or Macintosh computers. *The contents of this catalog are subject to change without notice.

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our Web site at www.shimadzu.com



SHIMADZU CORPORATION. International Marketing Division 3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan Phone: 81(3)3219-5641 Fax. 81(3)3219-5710 URL http://www.shimadzu.com

Printed in Japan 4199-02801-50AIT

